



GaussViewW Demo Installation Instructions

Thank you for choosing to review *GaussViewW*, a graphical user interface designed to work with *Gaussian 09*. Since you have received a demo version, you will need to contact Gaussian, Inc. to get a time-limited key for your trial period. When you are ready to begin the trial period, please send an email requesting a key to custserv@gaussian.com. The message must contain your CD serial number. The key and instructions for entering it will be sent by return email. Once you have completed your examination of the product, please contact us to arrange for purchase:

Email: orders@gaussian.com
Telephone: 203-284-2501
Fax: 203-284-2521

Installing the Demo Software

The distribution CD in this package is designed to install *GaussViewW* automatically on supported Windows systems. Simply insert the CD, and follow the installation directions on the screen. If the installation does not start automatically, you may start it manually selecting **Start=>Run** and then entering **D:\SETUP.EXE** into the resulting dialog (substitute a different drive letter if your CD-ROM drive is not D:).

The default installation directory for GVW is the same as your *Gaussian 09W* distribution, typically *C:\G09*. Using the same directory for both products simplifies GVW's interactions with the *Gaussian 09W* product.

During the installation you will be asked to supply a serial number, which you will find on the enclosed registration card. If your site has purchased a site license for *GVW*, your serial number is coded for multiple installations, and so the same serial number will work on each machine.

Note that although *GVW* will process *G03W* output, it will work best with *G09W*. If you have *G03W*, contact Gaussian, Inc. for information about upgrading to *G09W*.